





BIO-MEDICARE

BIOMEX 1010 FR-D-SA/RA

100 MA FULLWAVE RECTIFIED DIAGNOSTIC X-RAY SYSTEM




BIOMEX 1010 FR-D UNIT SA/RA



BASIC RADIOGRAPHY SYSTEM



BIO-MOBILE



MOTOREX TABLE

BIOMEX 1010 FR-D-SA/RA

TWO PULSE SINGLE TANK DIAGNOSTIC X-RAY UNIT

APPLICATION

BIOMEX 1010 FRD X-RAY UNIT SITUABLE FOR RADIOGRAPHIC AND ROUTINE FLUROSCOPY WORK

FEATURES

FULL WAVE RECTIFICATION

Full wave rectification Provides Maximum X-ray Output from available Power sources

LINE VOLTAGE COMPOSITION

Line voltage composition in 8 steps from 200 to 270 voltage

ELECTRONIC TIMER

Electronic timer 0.02 to 5 Sec. in 23 steps

DIGITAL DISPLAY OF MA-KV & MA-S

KV output major KV 40 KV to 90 KV in 10 KV step

KV output minor KV 0 KV to 10 KV in 2 KV step

MA output 40 MA to 100 MA in 4 step

Column Stand

Floor to Ceiling counter balanced column stand for X-ray tube having movements in all directions with locks. 360°. column rotation. "U" arm attachment & straight arm attachment

5 Position Table

5 position manual table with counter balanced motorised bucky with high ratio grid & stainless steel cassette tray take cassette upto 14" X 17" size. Horizontal bucky table with grid.

Table Dig. 15 trendenhenbug , 0°,30°,60°,90° manually.





BIO-MEX 1010 FR-D

SPECIFICATION	BIO-MEX1010 SA	BIO-MEX1010 RA	BIO-MEX 6010 FR-D
KV. OUTPUT	40 to 100 KV in 2 KV steps	40 to 100 KV in 2 KV steps	40 to 100 KV in 2 KV steps
MA. OUTPUT	40 to 100 MA in 4 steps	40 to 100 MA in 4 steps	20 to 60 MA in 4 steps
TIMER	0.02 to 5 Sec. in 23 steps	0.02 to 5 Sec. in 23 steps	0.04 to 5 Sec. in 23 steps
POWER INPUT	440 Volts 3 Phase (10KVA)	440 Volts 3 Phase (10KVA)	200 to 270 Volts Single Phase (15 Amp.)
X-RAY TUBE	BEL - DSA - 3 or Imported	BEL - DRA - 6 or Imported	BEL - DSA - 3 or Imported

SPOT FILM DEVICE

• Sturdy Spot film device equipped with 15" X 15" fine line X-ray grid. • Single inbuilt coning mechanism for easy operation • Magnetic locks • Cough guard & stay radiation guard • Accepts all standard size cassettes



SPOT FILM DEVICE

RADIATION PROTECTION CONFIRMED TO AERB STANDARD

NOTE: DUE TO OUR CONTINUOUS PRODUCT IMPROVEMENT PROGRAM, SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE

MANUFACTURED BY:

BIO MEDICARE IMAGING CO.

F-18 & 19, UDYAG VIHAR, INDUSTRIAL ESTATE,
VITHALWADI (WEST) (DIST.) THANE. PIN - 421 003. (MUMBAI)
TEL. 0251(911) - 543332. TELE FAX. 0251(911) - 543344.

MARKETED BY: