





BIO-MEDICARE

MICRO PROCESSOR BASED DIGITAL PORTABLE X-RAY UNITS




BIOMEX 1570-D

**Ultra Weight Portable Machine
Different Three Models**




BIOMEX 3085 FR-D

**Digital Micro Processor
Technology**



BIORATCHE STAND



BIOMEX 5085 FR-D

STANDARD FEATURES

- Truly Portable Carry away box type Machine
- Voltage Over Load Protection
- Tape built in FFD Measurement
- Micro Processor Based Digital Timer with Accuracy of 0.01 Sec.
- Exposure Time Can be set by soft touch buttons
- Full Proof Radiation Protection
- Line Voltage Compensation
- Hand Operated switch is provided for X-Ray exposure in single touch

• BIOMEX 1570 - D

Ultra Light Weight and Radiography Examination Done on Domestic Power Point

• BIOMEX 3085 FR-D

Variable KV/MA for Fine Radiography

• BIOMEX 5085 FR - D

Highest Capability Through Full Wave Rectification

SPECIFICATIONS :	BIOMEX 1570 - D	BIOMEX 3085 FR - D	BIOMEX 5085 FR - D
OUTPUT	15MA at 70 KV	30 MA at 50 KV 20 MA at 68 KV 15 MA at 85 KV	50 MA at 50 KV 40MA at 65 KV 35 MA at 65 KV 30 MA at 85 KV
TIMER	MICROPROCESSOR 0.1 to 6 Sec. 60 steps	MICROPROCESSOR 0.1 to 6 Sec. 60 steps	MICROPROCESSOR 0.1 to 6 Sec. 60 steps
TUBE	0.8 MM FOCUS TOSHIBA or Equivalent	1.8 MM FOCUS BEL DSA - 1 or Equivalent	1.8 MM FOCUS BEL DSA - 1 or Equivalent
INPUT COMPENSATION	230 V. 5 AMP. 210 TO 250 V.	230 V. 15 AMP. 210 TO 250 V.	230 V. 210 TO 250 V.
BEAM LIMITATION	Cone Centering Device	Cone Centering Device	Light Beam Collimator
WEIGHT	10 Kgs.	14 Kgs.	17 Kgs.

STANDARD ACCESSORIES :

Power Supply Cable Exposure Switch with Cable Carrying Case Beam Limitation Device

STANDS :

Bio-Light : Light Weight (12 Kgs.) Dismantable Stand

Bio Rathe - Sturdy Convenient Rathe Typestand

RADIATION PROTECTION CONFIRMED TO AERB STANDARD

NOTE: DUE TO OUR CONTINUOUS PRODUCT IMPROVEMENT PROGRAM, SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE

Manufactured by :

BIOMICARE IMAGING CO.

F-18, Udyog Vihar Ind. Estate, Vithalwadi (West)

Dist. Thane. Pin-421003 (Mumbai)

Tel. 95-251- (0251) 543332

Telefax : 95-251- (0251) 543344

MARKETED BY: